




Straightshot® M4 Microdebrider Procedure Guide



BLADES AND BURS FOR NASAL, SINUS, LARYNGEAL, AND BRONCHIAL SURGERY

1997-PRESENT

FOUR GENERATIONS OF POWER CONSOLES AND HANDPIECES



XPS®



XPS® 2000



XPS® 3000

Did you know?

Powered sinus debridors were developed and patented by Medtronic Xomed in 1997.

Powered surgery is firmly established for successfully managing both simple and complex ENT disorders. Since 1997, Medtronic ENT has developed a comprehensive range of powered surgery tools for performing minimally invasive ENT procedures in adults and children with greater precision and speed.

The current generation of powered ENT tools includes the IPC® Integrated Power Console and Straightshot® M4 Microdebrider. This dynamic combination provides:

- The widest range of application-specific blades and burs
- Blade tips that rotate 360°
- Factory-calibrated blades for surgical navigation
- IntelliFlow™ remote control irrigation
- Seamless integration with other Medtronic ENT devices




Straightshot® Magnum®



IPC® Integrated Power Console




Straightshot® M4 Microdebrider with
Rotatable Blades




POLYP REMOVAL

The advantages of the Straightshot® M4 Microdebrider include:¹

- **Excellent control and precision for resection of soft tissue**
- **Preservation of surrounding tissue because only tissue suctioned into the opening is cut and removed**
- **Enhanced visibility due to the suction removal of blood and debris from the surgical field**




Tricut® Blade - 1884004



	Product No.	Description	Comments
	1884004HR	Tricut Sinus Blade, M4 Rotatable, 4.0 mm – 5/box	3,000 - 5,000 rpm – oscillation
	1883504HR	Tricut Sinus Blade, M4 Rotatable, 3.5 mm – 5/box	3,000 - 5,000 rpm – oscillation

TURBINOPLASTY


Clinical research has shown significantly better long-term results with powered inferior turbinoplasty.²⁻⁷ The innovative Inferior Turbinate (IT) Blade has a patented elevator tip to facilitate insertion into the inferior turbinate and the creation of a submucosal pocket for bulk reduction. The oscillating IT blade is designed to quickly and easily debulk hypertrophied tissue within the turbinate.



Better than Coblation®

Visual Analog Scale (VAS) Scores

10 - The most severe symptoms



Key

Radiofrequency-Assisted Inferior Turbinoplasty -----
 Microdebrider-Assisted Inferior Turbinoplasty -----
 On the VAS Scale: 0 = No symptoms 10 = The most severe symptoms

Visual Analog Scale (VAS)

A subjective patient questionnaire that evaluates their perception of his or her health; in this case, pertaining to nasal obstruction, sneezing, rhinorrhea, and snoring. Answers usually range from zero (no symptoms) to 10 (the most severe symptoms).²


In the first three-year study directly comparing long-term efficacy of the two techniques, patients treated with Medtronic's microdebrider experienced significantly better long-term outcomes than those who received treatment with Coblation® radiofrequency.²

When results in the two groups were compared, none of the values for the assessed variables were significantly different at 6 months. At 1, 2, and 3 years, however, all results for both the subjective and objective variables were significantly better in the microdebrider group. Neither group experienced significant complications from surgery.

Product No.	Description	Comments
1882040HR	Inferior Turbinate Blade, M4 Rotatable, 2.0 mm – 5/box	60 - 3,000 rpm - oscillation
1882940HR	Inferior Turbinate Blade, M4 Rotatable, 2.9 mm – 5/box	60 - 3,000 rpm - oscillation



ETHMOIDECTOMY AND SPHENOIDOTOMY



One of the advantages of the Straightshot® M4 Microdebrider is being able to **use both straight and curved blades in a single procedure.**

Whether it's an anterior or total ethmoidectomy to manage chronic sinusitis or sinonasal polyps, the versatility of the microdebrider makes it possible to resect delicate ethmoidal tissue and bone to gain access to the sphenoid sinus.

The Straightshot M4 Microdebrider is frequently used in sinonasal polyposis that recur following an ethmoidectomy because it's associated with less bleeding after removal of the regrown polyps.¹



M4 Blades with Automated EM Tracking

The new M4 Blades with Automated EM Tracking deliver unparalleled convenience and technology integration. They're the first and only blades that are factory-calibrated for navigation, right out of the box. Attach the blade to the M4 Microdebrider and the Fusion™ system, and start navigating.



12°




1884012HR

40°





1884006EM

	Product No.	Description		Comments
360°	1884012HR	RAD® 12 Curved Sinus Blade, M4 Rotatable, 4.0 mm – 5/box	12°	3,000 - 5,000 rpm – oscillation
360°	1884006HR	RAD 40 Curved Sinus Blade, M4 Rotatable, 4.0 mm – 5/box	40°	3,000 - 5,000 rpm – oscillation
360°	1884016HR	RAD 60 Curved Sinus Blade, M4 Rotatable, 4.0 mm – 5/box	60°	3,000 - 5,000 rpm – oscillation
360°	1884012EM	RAD 12 Curved Rotatable Blade with Automated EM Tracking, 4.0 mm – 1/ea.	12°	3,000 - 5,000 rpm – oscillation
360°	1884006EM	RAD 40 Curved Rotatable Blade with Automated EM Tracking, 4.0 mm – 1/ea.	40°	3,000 - 5,000 rpm – oscillation
360°	1884080EM	Tricut® 13.0 cm Straight Rotatable Blade with Automated EM Tracking, 4.0 mm – 1/ea.		3,000 - 5,000 rpm – oscillation



Uncinectomy

UNCINECTOMY AND MAXILLARY ANTROSTOMY




Maxillary antrostomies are made easier with the M4 rotatable blades

With the Straightshot® M4's straight or curved blades, you can perform all of the following procedures without changing equipment:

- **Delicately resect the uncinate process**
- **Open the maxillary ostium to enable visual inspection of the maxillary cavity**
- **Curved blades (40, 60, 90, and 120 degree) enable resection of polyps and/or disease as necessary**

Precise resection of mucosa facilitates better healing of mucosal edges and ensures long-lasting patency.¹



	Product No.	Description		Comments
	1884012HR	RAD® 12 Curved Sinus Blade, M4 Rotatable, 4.0 mm – 5/box		3,000 - 5,000 rpm – oscillation
	1884006HR	RAD 40 Curved Sinus Blade, M4 Rotatable, 4.0 mm – 5/box		3,000 - 5,000 rpm – oscillation
	1884016HR	RAD 60 Curved Sinus Blade, M4 Rotatable, 4.0 mm – 5/box		3,000 - 5,000 rpm – oscillation
	1883519HR	RAD 90 Curved Sinus Blade, M4 Rotatable, 3.5 mm – 3/box		2,000 - 3,000 rpm – oscillation
	1883517	RAD 120 Curved Sinus Blade, 3.5 mm – 3/box		1,500 - 3,000 rpm – oscillation

TREPHINATION OF THE FRONTAL SINUS

Mini-Trephination Set





Designed for use with the Straightshot® M4 Microdebrider, the Mini-Trephination Set allows trephination through the anterior face of the frontal sinus for patients who require frontal sinus irrigation.

Product No.	Description	Comments
1892000*	Mini-Trephination Set – 1/box	6,000 rpm – right rotation

Set Includes:

2.0 mm Drill (1882900), Drill Guide with Skin Protector (1892001), Guide Pin (1892002), Irrigation Cannula (1892003), Instrument Tray (3717005)

**This set is only available in the USA. For customers outside the USA, please order the above items individually.*



DACRYOCYSTORHINOSTOMY

1882569HS

20°


15°

1884068HS



The Straightshot® M4 Microdebrider and DCR Burs are well-suited for performing endoscopic dacryocystorhinostomy (DCR). Simultaneous irrigation and aspiration of bone debris helps maximize visibility during drilling.

Product No.	Description	Comments
1884068HS	Curved DCR Bur 15°, 4.0 mm – 3/box	Up to 12,000 rpm
1882569HS	Curved Diamond DCR Bur 20°, 2.5 mm – 3/box	Up to 12,000 rpm




TONSILS AND ADENOIDS

Tonsil and Adenoid Blades for the PITA™ Technique (*Powered Intracapsular Tonsillectomy and Adenoidectomy*)

Clinical trials show that the PITA technique is safe, effective, and offers most patients less postoperative pain, a faster recovery, and a quicker return to normal diet.⁸⁻¹⁶

The T&A Set consists of two interchangeable blades, one curved (12°) and the other open towards the outside (40°).



1884008TA

Product No.	Description	Use	Comments
1884507	RADenoid® Adult Blade, 4.5 mm, 45° curved – 5/box	Adenoids in adults	1,500 rpm – oscillation
1884008	RADenoid Pediatric Blade, 4.0 mm, 40° curved – 5/box	Adenoids in children	1,500 rpm – oscillation
1884008TA	Tonsil (12°, 4.0 mm) & Adenoid (40°, 4.0 mm) Blade Set – 5/box	T&A	1,500 rpm – oscillation
1884013	Curved Tonsillectomy Blade, 12°, 4.0 mm – 5/box	Tonsillectomy	1,500 rpm – oscillation



CHOANAL ATRESIA

The Straightshot® M4 Microdebrider enables effective tissue reduction in simple soft tissue or complex boney choanal atresia procedures.¹⁷⁻²⁰




1882960




1882905



1883673HS



Choanal atresia in the right nasal fossa



Creation of a new opening using a 2.9 mm bur

Product No.	Description	Use	Comments
1882905	Straight Conical Silver Bullet® Blade, 2.9 mm – 5/box	Pediatric	3,000 - 5,000 rpm - oscillation
1882960	Round Straight Cutting Bur, 2.9 mm – 5/box	Pediatric	5,000 rpm - oscillation
1883262HS	Round Straight Cutting Bur, 3.2 mm – 3/box	Pediatric	12,000 rpm – right rotation
1883673HS	Choanal Atresia Bur, 4.0 mm, 15° curved – 3/box	Pediatric	12,000 rpm – rotation

Septoplasty and Rhinoplasty



Specialty Burs

1883212HS

1884566

1883260

The FeatherTouch® Rasp's unique design converts rotation to reciprocating movement, which is especially useful in rhinoplasty. This enables effortless bone removal and an appealing nasal profile.




1992208



1992210

Product No.	Description	Use	Comments
1883212HS	Septoplasty Bur, 3.2 mm, 12° Curved – 3/box	Septoplasty	12,000 rpm – right rotation
1884566	RhinoBur®, 4.5 mm – 3/box	Rhinoplasty	4,000 – 6,000 rpm – right rotation
1883260	Tardy MicroBur®, 3.2 mm – 5/box	Rhinoplasty	3,000 - 5,000 rpm – right rotation
1922005	FeatherTouch® Converter – 1/box	Rhinoplasty - hump reduction	Together with the above-mentioned rasps
1992208	FeatherTouch Rasp, Coarse (width: 8.0 mm) – 2/box		3,000 - 5,000 rpm – mounted on the above-mentioned converter
1992210	FeatherTouch Rasp, Fine (width: 8.0 mm) – 2/box		




Larynx and Airways

Specialty Blades

- Papillomatosis, subglottic, and tracheal stenosis
- Polyps, edema, or tumors in the airways

With the Straightshot® M4 Microdebrider, tissue can be removed quickly and precisely to restore the airway as much as possible without compromising the larynx. The Skimmer® blade enables **safe, fast, and precise removal of airway lesions** in recurrent respiratory papillomatosis (RRP), eliminating thermal burn and other risks of laser techniques.²¹

Further indications may include unblocking of the larynx and bronchi and phonosurgery for polypoid degeneration.



Product No.	Description	Use	Comments
1882925	Skimmer® Angle-Tip Blade, 15°, 2.9 mm, L: 18.5 cm – 3/box	Larynx Papilloma	60 – 500 rpm - oscillation
1882925HRE	Skimmer Angle-Tip Blade, 15°, M4 Rotatable, 2.9 mm, L: 18.0 cm – 1/box	Larynx Papilloma	60 – 500 rpm - oscillation
1884023	Skimmer Angle-Tip Blade, 15°, 4.0 mm, L: 22.5 cm – 3/box	Larynx Papilloma	60 – 500 rpm - oscillation
1882923HRE	Skimmer Angle-Tip Blade, 15°, M4 Rotatable, 2.9 mm, L: 22.0 cm – 1/box	Larynx Papilloma	60 – 500 rpm - oscillation
1884024	Skimmer Angle-Tip Blade, 15°, 4.0 mm, L: 27.5 cm – 3/box	Larynx Papilloma	60 – 500 rpm - oscillation
1882924HRE	Skimmer Angle-Tip Blade, 15°, 2.9 mm, L: 27.0 cm – 1/box	Larynx Papilloma	60 – 500 rpm - oscillation
1884020	Straight Tricut® Laryngeal Blade, 4.0 mm, L: 22.5 cm – 3/box	Debulking	60 – 1,200 rpm - oscillation
1884030HRE	Tricut Angle-Tip Blade, 15°, M4 Rotatable, 4.0 mm, L: 22.0 cm – 1/box	Larynx Debulking	60 – 1,200 rpm - oscillation
1884031HRE	Tricut Angle-Tip Blade, 18°, M4 Rotatable, 4.0 mm, L: 27.0 cm – 1/box	Trachea Debulking	60 – 1,200 rpm - oscillation
1884033HRE	Curved Tracheal Blade, M4 Rotatable, 15°, 4.0 mm, L: 37.0 cm – 1/box	Trachea and bronchi	60 – 1,200 rpm – oscillation
1884035HRE	Curved Bronchial Blade, M4 Rotatable 4.0 mm, L: 45.0 cm – 1/box	Bronchi	60 – 2,000 rpm - oscillation

References

1. Gendeh BS, Long Y-T, Misiran K. Antrochoanal polyps: clinical presentation and the role of powered endoscopic polypectomy. *Asian J Surg* 2004; 27: 22-5.
2. Liu C-M, Tan C-D, Lee F-P, Lin K-N, Huang H-M. Microdebrider-assisted versus radiofrequency-assisted inferior turbinoplasty. *Laryngoscope* 2009; 119: 414-8.
3. Yañez C. Inferior turbinate debridement technique: Ten-year results. *Arch Otolaryngol Head Neck Surg* 2008; 138:170-175.
4. Berger G, Ophir D, Pitro K, Landsberg R. Histopathological changes after Coblation® inferior turbinate reduction. *Arch Otolaryngol Head Neck Surg* 2008;134(8):819-823.
5. Huang TW, Cheng P. Changes in nasal resistance and quality of life after endoscopic microdebrider-assisted inferior turbinoplasty in patients with perennial allergic rhinitis. *Arch Otolaryngol Head Neck Surg* 2006; 132: 990-993.
6. Sacks R, Thornton MA, Boustred RN. Modified endoscopic turbinoplasty—long-term results compared to submucosal electrocautery and submucosal powered turbinoplasty. Presented at: American Rhinologic Society Spring Meeting; May 13–16, 2005; Boca Raton, FL.
7. Atef A, Mosleh M, El Bosraty H, El Fatah GA, Fathi A. Bipolar radiofrequency volumetric tissue reduction of inferior turbinate: Does the number of treatment sessions influence the final outcome? *Am J Rhinol* 2006; 20: 25-31.
8. Koltai PJ, Kalathia AS, Stanislaw P, and Heras HA. Power-Assisted Adenoidectomy. *Archives of Otolaryngology Head & Neck Surgery*, July 1997; 123: 685-88.
9. Stanislaw P, Koltai PJ, and Feustel PJ. Comparison of Power-Assisted Adenoidectomy vs. Adenoid Curette Adenoidectomy. *Archives of Otolaryngology Head & Neck Surgery*, July 2000; 126: 845-49.
10. Murray, LN, Fitzpatrick P, and Guarisco JL. Powered Partial Adenoidectomy: A Clinical Trial. *Archives of Otolaryngology Head & Neck Surgery*, July 2002; 128: 792-796.
11. Rodriguez K, Murray JN, and Guarisco JL. Powered-Assisted Partial Adenoidectomy. *Laryngoscope*, August 2002; 112: 26-28.
12. Koltai PJ, Chan J, and Younes A. Powered-Assisted Adenoidectomy: Total and Partial Resection. *Laryngoscope*, August 2002; 112: 29-31.
13. Havas, Thomas, MD and Lowinger, David, MD. Obstructive Adenoid Tissue: An Indication for Powered-Shaver Adenoidectomy. *Archives of Otolaryngology Head & Neck Surgery*, July 2002; 128: 789-791.
14. Koltai PJ, Solares A, Mascha EJ, and Xu M. Intracapsular Partial Tonsillectomy for Tonsillar Hypertrophy in Children. *Laryngoscope*, August 2002; 112: 17-19.
15. April M, Ward R, Bent J. Power-Assisted Adenoidectomy in the Treatment of Chronic Otitis Media with Effusion. Poster Presentation at American Society of Pediatric Otolaryngology, May 4, 2003, Nashville, TN.
16. Koltai PJ, Solares CA, Abelson TI, Hirose K, Krakovitz PR, Reilly JS, Cook SP, April MM, Ward RF, Bent JP, Bhayani R, Reilly B, Weiss G. The Safety and Efficacy of Intracapsular (Partial) Tonsillectomy (IT) for Pediatric Obstructive Sleep Disordered Breathing (OSDB). Presented at the American Society of Pediatric Otolaryngology (ASPO), May 5, 2003, Nashville TN.
17. Froehlich P, Pransky SM, Fontaine P, Stearns G, Morgon A. Minimal endoscopic approach to subperiosteal orbital abscess *Arch Otolaryngol Head Neck Surg* 1998 Jun; 124(6):711.
18. McLeod IK, Brooks DB, Mair EA. Revision choanal atresia repair. *Int J Pediatr Otorhinolaryngol* 2003; 67: 517-24.
19. Rombaux P, de Toeuf C, Hamoir C, Eloy P, Bertrand B, Veykemens F. Transnasal repair of unilateral choanal atresia. *Rhinology* 2003; 41: 31-6.
20. Schoem SR. Transnasal endoscopic repair of choanal atresia: why stent? *Otolaryngol Head Neck Surg* 2004; 131: 362-6.
21. El-Bitar M, Zalzal GH. Powered instrumentation in the treatment of recurrent respiratory papillomatosis: an alternative to the carbon dioxide laser. *Arch Otolaryngol Head Neck Surg* 2002;128:425-8.

For further information, please call Medtronic ENT at 800.874.5797 or 904.296.9600. You may also consult our website at www.MedtronicENT.com.

Medtronic ENT

Medtronic USA, Inc.

6743 Southpoint Drive N
Jacksonville, FL 32216
USA
www.MedtronicENT.com
Toll free: (800) 874-5797
Fax: (800) 678-3995

International Telephone Numbers

Adriatic Reg. Office 385-1-488-1120	Italy 39-02-24137-324
Australia 1800-668-670	Japan 81-6-4795-1506
Baltic Reg. Office 37-1-67560226	Korea 82-2-3404-3600
Belgium 32-2456-09-09	Lebanon 961-1-370-670
Canada 1800-217-1617	Luxembourg 32-2456-09-09
China 86-21-50800998	Netherlands 31-45-566-8800
Czech Republic 420-2-9657-9580	Poland 48-22-465-6942
France 33-470-679-800	Russian Fed. 7-495-580-73 77
Germany 49-2159-8149-209	Singapore 65-6776-6255
Greece 30-210-67-79-099	South Africa 27-11-466-1820
Hong Kong 852-2919-1312	Spain 34-91-625-05-40
Hungary 36-30-5052987	Taiwan 886-2-2183-6000
India 91-22-26836733	UK 44-1923-205-166
Israel 972-9-972-4400	USA 1-904-296-9600



Medtronic