


Foot Switching Bipolar Forceps, Non-Sterile, Reusable

Steam/EtO sterilizable

Only sold in selected countries




Tip shown actual size

E4051-CT

1 EACH

Scoville-Greenwood Bayonet Forceps†

- Insulated
- Smooth tips
- Total length: 19,7 cm (7,75 in.)
- Tip: 1,5 mm



Tip shown actual size

E4054-CT

1 EACH

Cushing Bayonet Forceps†

- Insulated
- Smooth tips
- Total length: 19,1 cm (7,5 in.)
- Tip: 2,0 mm



Tip shown actual size

E4052-CT

1 EACH

Cushing Forceps†

- Insulated
- Smooth tips
- Total length: 15,2 cm (6 in.)
- Tip: 0,7 mm



Tip shown actual size

E4055-CT

1 EACH

Jewelers Forceps, Straight Tip†

- Insulated
- Smooth tips
- Total length: 10,2 cm (4 in.)
- Tip: 0,4 mm




Tip shown actual size

E4053-CT

1 EACH

Cushing Forceps†

- Insulated
- Smooth tips
- Total length: 17,8 cm (7 in.)
- Tip: 1,5 mm



Tip shown actual size

E4057-CT

1 EACH

Gerald Bayonet Forceps†

- Insulated
- Smooth tips
- Total length: 19,1 cm (7,5 in.)
- Tip: 0,7 mm

† Requires E0512 single use cords or QI2401 or E0020V reusable cords

Foot Switching Bipolar Forceps, Non-Sterile, Reusable

Steam/EtO sterilizable

Only sold in selected countries



Tip shown actual size

E4058-CT

1 EACH

Cushing Bayonet Forceps

- Insulated
- Smooth tips
- Total length: 19,1 cm (7,5 in.)
- Tip: 0,7 mm




Tip shown actual size

E4062-CT

1 EACH

Curved Iris Forceps

- Insulated
- Smooth tips
- Total length: 8,9 cm (3,5 in.)
- Tip: 0,5 mm



Tip shown actual size

E4059-CT

1 EACH

Hardy Bayonet Forceps with Stops*

- Insulated
- Smooth tips
- Total length: 20,9 cm (8,25 in.)
- Tip: 0,5 mm



Tip shown actual size

E4073-CT

1 EACH

Titanium Bayonet Forceps

- Insulated
- Smooth tips
- Total length: 22,2 cm (8,75 in.)
- Tip: 1,0 mm



Tip shown actual size

E4060-CT

1 EACH

Semkin Forceps with Stops*

- Insulated
- Smooth tips
- Total length: 14 cm (5,5 in.)
- Tip: 0,5 mm

*Stops prevent scissoring of the tip

Bipolar Coagulation Forceps, Non-Sterile, Reusable



E700150

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Straight
- Pointed



E700157

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Angled
- Tip: 1 mm



E700151

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Straight
- Tip: 1 mm



E700158

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Short angled
- Tip: 0,7 mm



E700155

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Angled
- Pointed
- Delicate tines



E700167

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Bayonet
- Tip: 1 mm

Bipolar Coagulation Forceps, Non-Sterile, Reusable



E700170 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Straight
- Pointed



E700171 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Straight
- Tip: 1 mm



E700175 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Angled
- Tip: 1 mm

E700176 1 EACH

Bipolar Coagulation Forceps

Same specifications except:

- Tip: 2 mm



E700181 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Bayonet
- Tip: 1 mm



E700184 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Bayonet
- Angled downwards
- Tip: 1 mm



E700185 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Bayonet
- Angled upwards
- Tip: 1 mm

Bipolar Coagulation Forceps, Non-Sterile, Reusable



E700191 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 23 cm (9,1 in.)
- Bayonet
- Tip: 0,7 mm

E700192 1 EACH

Bipolar Coagulation Forceps

Same specifications except:

- Tip: 1,2 mm



E700195 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 23 cm (9,1 in.)
- Angled
- Pointed
- Delicate tines



E700240 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 10,5 cm (4,25 in.)
- Straight
- Pointed



E700241 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 10,5 cm (4,25 in.)
- Straight
- Blunt



E700243 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 10,5 cm (4,25 in.)
- Straight
- Blunt
- Tip: 1 mm



E700244 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 10,5 cm (4,25 in.)
- Short angled
- Pointed

Bipolar Coagulation Forceps, Non-Sterile, Reusable



E700204 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 14,5 cm (5,75 in.)
- Tip: 1 mm



E700214 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 14,5 cm (5,75 in.)
- Tip: 2 mm



E700206 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Tip: 1 mm



E700216 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Tip: 2 mm



E700209 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Tip: 1 mm



E700219 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Tip: 2 mm



E700212 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 22 cm (8,75 in.)
- Tip: 1 mm



E700222 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 22 cm (8,75 in.)
- Tip: 2 mm

Bipolar Coagulation Forceps, Non-Sterile, Reusable



E700245

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 10,5 cm (4,25 in.)
- Long angled
- Pointed



E700284

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Bayonet
- Angled downwards
- Tip: 0,7 mm



E700246

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 10,5 cm (4,25 in.)
- Long angled
- Blunt



E700292

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 25 cm (9,8 in.)
- Bayonet
- Tip: 1,2 mm



E700247

1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 11 cm (4,33 in.)
- Long angled
- Blunt
- Tip: 1 mm

Bipolar Coagulation Forceps, Non-Sterile, Reusable



E700331 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 11 cm (4,5 in.)
- Straight
- Tip: 1,2 mm



E700352 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 22 cm (8,75 in.)
- Straight
- Tip: 2,2 mm



E700346 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Straight
- Tip: 2,2 mm



E700355 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 25 cm (10 in.)
- Straight
- Tip: 2,2 mm



E700349 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Straight
- Tip: 2,2 mm



E700357 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 30 cm (12 in.)
- Straight
- Tip: 2,2 mm

Bipolar Coagulation Forceps, Non-Sterile, Reusable



E700376 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Angled
- Tip: 2,2 mm



E700760¹ 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 16,5 cm (6,5 in.)
- Bayonet
- Tips: 1,1 mm



E700382 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 19 cm (7,5 in.)
- Bayonet
- Tip: 2,2 mm



E700764¹ 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 20 cm (7,9 in.)
- Bayonet
- Tip: 1,1 mm



E700768¹ 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 23 cm (9,1 in.)
- Bayonet
- Tip: 1,1 mm



E700770¹ 1 EACH

Bipolar Coagulation Forceps

- Use with E360XXX reusable cord
- Total length: 25 cm (9,8 in.)
- Bayonet
- Tip: 1,1 mm

1. Made to order