


ENT INTERNATIONAL PRODUCT CATALOG




HYALURONIC ACID-BASED DRESSING, PACKING, AND STENTING PRODUCTS

Medtronic

MeroGel™ Nasal Dressing/Sinus Stent/Otologic Dressing or Packing

- A packing material made from HYAFF™, an esterified form of hyaluronic acid


Product	Description	Qty
1517100	MeroGel Otologic Packing, 1 cm x 5 cm	1
1517000	MeroGel Nasal Dressing/Sinus Stent and Otologic Packing, 4 cm x 4 cm	1
1517002	MeroGel Nasal Dressing/Sinus Stent and Otologic Packing, 4 cm x 4 cm	2
1518000	MeroGel Injectable Bioresorbable Nasal Dressing/Sinus Stent and Otologic Gel (includes one syringe and one flexible catheter)	1
1520000	MeroPack Bioresorbable Nasal Dressing and Sinus Stent	2
1520001	MeroPack Bioresorbable Nasal Dressing and Sinus Stent	6



EpiDisc™ Otologic Lamina

- An otologic device made from HYAFF™, an esterified form of hyaluronic acid
- Pre-cut, 8 mm spherical discs of EpiFilm™ material

Product	Description	Qty
1417100	EpiDisc Otologic Lamina, 8 mm Disc	2



EpiDisc™ Tympanic Membrane Perforation Patch Kit

- Consists of one 8 mm EpiDisc matrix device and a 2 cm x 2 cm square of MeroGel™ QuickGel™ Otologic Dressing, a fast-hydrating form of MeroGel


Product	Description	Qty
1417151	EpiDisc Tympanic Membrane Perforation Patch Kit	1



EpiFilm™ Otologic Lamina

- An otologic device made from HYAFF™, an esterified form of hyaluronic acid

Product	Description	Qty
1417000	EpiFilm Otologic Lamina, 2.5 cm x 2.5 cm	1




HYALURONIC ACID-BASED DRESSING, PACKING, AND STENTING PRODUCTS

POST-OP PACKING, SPLINTS, AND EPISTAXIS TREATMENT PRODUCTS

Medtronic

Merocel Standard Nasal Packing

Product	Description	Length	Width	Height	Qty
400400	Merocel without String	4.5 cm	1.5 cm	2 cm	20
440400	Merocel with String	4.5 cm	1.5 cm	2 cm	10
450400	Merocel hemoX	4.5 cm	1.5 cm	2 cm	10
470401	Merocel 2000 with String	4.5 cm	1.5 cm	2 cm	10
400402	Merocel without String	8.0 cm	1.5 cm	2 cm	20
440402	Merocel with String	8.0 cm	1.5 cm	2 cm	10
450402	Merocel hemoX	8.0 cm	1.5 cm	2 cm	10
400403	Merocel without String, Single Compression	8.0 cm	1.5 cm	2 cm	20
470404	Merocel 2000 with String	8.0 cm	1.5 cm	2 cm	10
400405	Merocel without String, Single Compression	8.0 cm	1.5 cm	2 cm	20
470405	Merocel 2000 with String, with Airway	8.0 cm	1.5 cm	2 cm	10
470408	Merocel 2000 with String, with Airway	4.5 cm	1.5 cm	2 cm	10
400409	Merocel without String, with Airway	4.5 cm	1.5 cm	2 cm	20
440409	Merocel with String, with Airway	4.5 cm	1.5 cm	2 cm	10
450409	Merocel hemoX with Airway	4.5 cm	1.5 cm	2 cm	10
400410	Merocel without String, with Airway	8.0 cm	1.5 cm	2 cm	20
440410	Merocel with String, with Airway	8.0 cm	1.5 cm	2 cm	10
450410	Merocel hemoX with Airway	8.0 cm	1.5 cm	2 cm	10





Merocel Epistaxis Nasal Packing

- Available with or without drawstring

Pope Flex-Pak™

- Epistaxis packing with pinholes and drawstring

Product	Description	Length	Width	Height	Qty
400406	Merocel Pope without String	10.0 cm	1.5 cm	2.5 cm	10
440406	Merocel Pope with String	10.0 cm	1.5 cm	2.5 cm	10
450406	Merocel hemoX Pope	10.0 cm	1.5 cm	2.5 cm	10
460406	Merocel Pope Flex-Pak™	10.0 cm	1.5 cm	2.5 cm	10
400416	Merocel Staxi-Stat™ Pack without String	-	-	-	10 pairs
440416	Merocel Staxi-Stat™ Pack with String	-	-	-	10 pairs
450416	Merocel hemoX Staxi-Stat™ with String	-	-	-	10
400421	Merocel Weimert with String	5.5 cm	1.5 cm	2.5 cm	20
450421	Merocel hemoX Weimert	5.5 cm	1.5 cm	2.5 cm	10



POST-OP PACKING, SPLINTS, AND EPISTAXIS TREATMENT PRODUCTS

All dimensions are approximate.

All dimensions are approximate.

Merocel Doyle Nasal Packing



Product	Description	Length	Thickness	Height	Qty
400411	Merocel without String, with Airway	8.0 cm	1.5 cm	3.0 cm	10
440411	Merocel with String, with Airway	8.0 cm	1.5 cm	3.0 cm	10
450411	Merocel hemoX with Airway	7.7 cm	1.5 cm	3.3 cm	10
400412	Merocel without String, Pediatric	5.5 cm	1.0 cm	2.5 cm	20
440412	Merocel with String, Pediatric	5.5 cm	1.0 cm	2.5 cm	10
400413	Merocel without String	8.0 cm	1.0 cm	3.0 cm	10
440413	Merocel with String	8.0 cm	1.0 cm	3.0 cm	10
450413	Merocel hemoX Slim-Profile	7.7 cm	1.0 cm	3.3 cm	10
400424	Merocel without String	8.0 cm	1.5 cm	3.0 cm	10
440424	Merocel with String	8.0 cm	1.5 cm	3.0 cm	10
450424	Merocel hemoX Regular	7.7 cm	1.5 cm	3.3 cm	10



Merocel Baron Sinus-Pak

- Three thicknesses



Product	Description	Length	Width	Height	Qty
400407	Merocel Baron Sinus-Pak with String, Regular	2.5 cm	1.2 cm	2.0 cm	10 pairs
450407	Merocel hemoX Baron Sinus-Pak Regular	2.5 cm	1.2 cm	2.0 cm	10 pairs
400427	Merocel Baron Sinus-Pak with String, Slim-Profile	2.5 cm	0.9 cm	2.0 cm	10 pairs
400531	Merocel Baron Sinus-Pak with String, Ultra-Slim Profile	2.5 cm	0.6 cm	2.0 cm	10 pairs
450531	Merocel hemoX Baron Sinus-Pak Ultra-Slim	2.5 cm	0.6 cm	2.0 cm	10 pairs
470531	Merocel 2000 Laminated Baron SinuStent™ Packing with String, Ultra-Slim Profile	2.5 cm	0.6 cm	2.0 cm	10 pairs



Merocel Kennedy Sinus-Pak

- Three thicknesses

Product	Description	Length	Width	Height	Qty
400422	Merocel Kennedy Sinus-Pak with String, Regular	3.5 cm	1.2 cm	1.2 cm	10 pairs
450422	Merocel hemoX Kennedy Sinus-Pak, Regular	3.5 cm	1.2 cm	1.2 cm	10 pairs
400426	Merocel Kennedy Sinus-Pak with String, Slim-Profile	3.5 cm	0.9 cm	1.2 cm	10 pairs
450426	Merocel hemoX Kennedy Sinus-Pak, Slim-Profile	3.5 cm	0.9 cm	1.2 cm	10 pairs
400530	Merocel Kennedy Sinus-Pak with String, Ultra-Slim Profile	3.5 cm	0.6 cm	2.0 cm	10 pairs
450530	Merocel hemoX Kennedy Sinus-Pak, Ultra-Slim	3.5 cm	0.6 cm	1.2 cm	10 pairs
470530	Merocel 2000 Laminated Kennedy SinuStent With String, Ultra-Slim Profile	3.5 cm	0.6 cm	1.2 cm	10 pairs




All dimensions are approximate.

All dimensions are approximate.

Merocel Toffel Sinus Packing

Product	Description	Length	Width	Height	Qty
400428	Merocel Toffel Slim Profile Sinus-Pak	4.5 cm	0.9 cm	2.0 cm	10 pairs
450428	Merocel hemoX Toffel Slim Profile Sinus-Pak	4.5 cm	0.9 cm	2.0 cm	10 pairs
400532	Merocel Toffel Ultra-Slim Profile Sinus-Pak	4.5 cm	0.6 cm	2.0 cm	10 pairs



Toffel Sinus-Pak Kit with Airway Tube

Product	Description	Length	Width	Height	Qty
400423	Toffel Sinus-Pak Kit with Airway Tube	—	—	—	10 kits


System Includes:

Sinus-Pak	4.5 cm	1.3 cm	2.0 cm
Nasal Packing	7.5 cm	1.5 cm	3.5 cm
Airway Tube	8.0 cm	—	—

Kennedy Intranasal Surgical Sponge (KISS) Intra-Op Throat Pack


Designed by David Kennedy, MD, Philadelphia, PA

Product	Description	Qty
400515	Kennedy Intranasal Surgical Sponge (KISS) Intra-Op Throat Pack	10



Oro-Pak Intra-Op Throat Pack

Product	Description	Qty
400520	Oro-Pak Intra-Op Throat Pack	10



All dimensions are approximate.

Merocel Pope Ear Wicks

Product	Description	Length	Width	Qty
400141	Merocel Pope	9.0 mm	15.0 mm	50
425141	Merocel Pope	9.0 mm	15.0 mm	25
400146	Merocel Pope	9.0 mm	24.0 mm	10
400144	Merocel Pope	12.0 mm	24.0 mm	10
450141	Merocel hemoX Pope	1.5 mm	0.9 mm	50



Merocel Pope Oto Wicks

Product	Description	Size	Qty
1425008	Merocel Pope Fenestrated	15.0 mm	50
1425010	Merocel Pope	24.0 mm	50



Merocel Schindler Ear Packing

Product	Description	Length	Thickness	Qty
400355	Merocel Schindler with Drawstring	12.0 mm	15.0 mm	10
400356	Merocel Schindler with Drawstring	9.0 mm	13.0 mm	10
400357	Merocel Schindler with Drawstring, Fenestrated	9.0 mm	13.0 mm	10
450355	Merocel hemoX Schindler	1.5 mm	1.2 mm	10
450356	Merocel hemoX Schindler Fenestrated	1.3 mm	0.9 mm	10



All dimensions are approximate.

Merocel Ambrus Ear Packing

Product	Description	Diameter	Length	Qty
400354	Merocel Ambrus	15 mm	20.00 mm	20

Merocel Ear Packs, Post-Op

Product	Description	Diameter	Length	Qty
400162	Merocel Ear Packs, Post-Op	12 mm	24.00 mm	10

Merocel Micro Spheres

Product	Description	Diameter	Thickness	Qty
400150	Merocel Micro Spheres	3 mm	1.00 mm	100
400151	Merocel Micro Spheres	4 mm	1.00 mm	100
400152	Merocel Micro Spheres	5 mm	1.25 mm	100



Medtronic ENT Epistat™ II Nasal Catheter


- Posterior balloon and Merocel™ Sponge combined
- Adjustable retaining ring
- Airway can be removed, leaving sponge in place

Product	Description	Length	Thickness	Height	Qty
400470	Catheter Only	5.5 cm	1.5 cm	2.5 cm	1

Medtronic ENT Epistat™ Nasal Catheter

- Two independently inflatable cuffs
- Hollow catheter
- Silicone construction
- Beveled posterior tip

Product	Description	Qty
1527031	Catheter Only	1



Tami Sinus Secretion Collector

- Sterile, disposable aspirator and specimen collector
- Developed in conjunction with Thomas A. Tami, MD

Product	Description	Qty
1514000	Tami Sinus Secretion Collector	10



All dimensions are approximate.

All dimensions are approximate.

Boomerang™ Turbinate Glove and Septal Splint

Product	Description	Qty
600500	Boomerang Turbinate Glove and Septal Splint	5 pairs

Oval Septal Splints

Product	Lenght	Thickness	Height	Qty
1523000	6,4 cm	1 mm	4,2 cm	5 pairs

Bridgemaster™ Nasal Splints

- Articulated nasal splint's flexible aluminum stabilizer adjusts inward
- Segmented foam pads

Product	Description	Qty
600302	Small Splint	5
600306	Large Splint	5

Doyle Bi-Valve Septal Splint

Product	Description	Qty
1524075	Doyle Bi-Valve Septal Splint	1 pair

Doyle II Intranasal Airway Splints

- Integrated airway
- Medical grade silicone
- Hole in anterior tip

Product	Description	Qty
1524050	Doyle II Intranasal Airway Splints	5 pairs
1524055	Doyle II Intranasal Airway Splints	1 pair



External Nasal Splints

- Lightweight aluminum lined with a thin layer of foam rubber

Product	Description	Qty
1528110	Small	25
1528111	Small	5
1528120	Medium	25
1528121	Medium	5
1528137	Large	25
1528138	Large	5



Goldsmith Inflation Airway Splint

- Inflatable cuff
- Integrated airway
- Medical grade silicone

Designed by Manning Goldsmith, III, MD, Savannah, GA

Product	Description	Qty
1524060	Goldsmith Inflation Airway Splint	1 pair



Goode Magne-Splint

- Cobalt-Samarium magnets ensure compression
- Medical grade silicone
- Pre-attached string

Designed by Richard L. Goode, MD, Palo Alto, CA

Product	Description	Qty
1516000	Regular	1 pair
1516010	Large	1 pair



All dimensions are approximate.


All dimensions are approximate.

Reuter Splint

- Constructed of fluoroplastic
- Pre-punched holes for easier suturing

Designed by S. Harold Reuter, MD, Houston, TX

Product	Description	Qty
1527010	Thin 0.25 mm	5 pairs
1527015	Thick 0.50 mm	5 pairs
1527020	Thick 0.50 mm, No Split	5 pairs
1527025	Thin 0.25 mm, No Split	5 pairs
1527000	Thin 0.25 mm	5 pairs
1527005	Thick 0.50 mm	5 pairs



Silicone Splint

- Septal splint folds
- Firm, yet flexible silicone
- Available to order individually or as kit
- Kit includes:
 - 2 Silicone Septal Splints (400561)
 - 2 Merocel™ Nasal Packs

Product	Description	Length	Thickness	Height	Qty
400561	Silicone Splint	75.00 mm	25.00 mm	1.00 mm	5 pairs
400562	Silicone Splint Kit				5 kits



Simple Splint


- Simple, inexpensive fluoroplastic splint
- Easily trimmed to fit

Product	Description	Length	Thickness	Height	Qty
1522000	Simple Splint	7.00 cm (2.75")	0.50 cm (0.02")	4.40 cm (1.75")	10 pairs



Thermasplints Kit

Product	Description	Length	Height	Qty
1529000	Small	1.25 cm	2.25 cm	10
1529010	Medium	1.53 cm	2.38 cm	10
1529020	Large	1.73 cm	2.63 cm	10



Nasal Septal Button

- Medical silicone construction for maximum compatibility

Product	Diameter	Qty
1524105	3 cm	1




Two-Part Adjustable Septal Button

- For non-surgical closure of nasal septal perforations
- Two-part button
- Medical silicone construction

Designed by Manning Goldsmith, III, MD, Savannah, GA

Product	Diameter	Qty
1524110	3 cm	1



All dimensions are approximate.

Instrument Wipes

- Removes fluid and debris from delicate instruments

Product	Description	Qty
2218525	3.50" x 3.50" (8.9 cm x 8.9 cm)	25
2218530	3.50" x 3.50" (8.9 cm x 8.9 cm), Adhesive-Backed	25

Medtronic ENT Silicone Sheeting

- Sheets are made of biocompatible silicone elastometer
- Sheets can be cut and trimmed as needed during procedures

Product	Description	Qty
1532500	5.00 cm x 5.00 cm x 0.25 mm (2" x 2" x .010")	10
1532505	5.00 cm x 5.00 cm x 0.13 mm (2" x 2" x .005")	10
1532510	5.00 cm x 5.00 cm x 0.50 mm (2" x 2" x .020")	10
1532520	5.00 cm x 5.00 cm x 1.00 mm (2" x 2" x .040")	10
1532530	5.00 cm x 5.00 cm x 1.60 mm (2" x 2" x .062")	10

Merocel Surgical Patties and Strips

- Made of soft, absorbent, fiber-free Merocel sponge
- 10 per packet, 10 packets per box



Product	Description	Qty
300501	1.3 cm x 1.3 cm (0.50" x 0.50")	100
300503	1.3 cm x 2.5 cm (0.50" x 1.00")	100
300505	1.3 cm x 5.0 cm (0.50" x 2.00")	100
300506	1.3 cm x 7.6 cm (0.50" x 3.00")	100
300507	2.5 cm x 2.5 cm (1.00" x 1.00")	100
300509	2.5 cm x 5.1 cm (1.00" x 2.00")	100
300510	2.5 cm x 7.6 cm (1.00" x 3.00")	100
300514	8.3 cm x 8.3 cm (3.25" x 3.25")	100
300515	0.7 cm x 7.6 cm (0.25" x 3.00")	100

Micro-Grid Silicone Metric Ruler

- Metric scale background
- Radiopaque medical grade silicone
- Sterile

Product	Description	Qty
4022800	Green, 30.0 mm x 15.0 mm x 0.4 mm with 1.0 mm Squares	10
4022805	Blue, 30.0 mm x 15.0 mm x 0.4 mm with 1.0 mm Squares	10

All dimensions are approximate.



ALAR™ Nasal Valve Stent


Designed in conjunction with Howard Stupak M.D.

Collapse of the upper lateral cartilage is associated with narrowing of the internal nasal valve, reducing nasal airflow.

The ALAR™ Nasal Valve Stent is intended for use in supporting the lateral nasal wall during the post-operative healing period following nasal valve surgery. This procedure may be performed independently or in conjunction with other procedures such as septoplasty and turbinoplasty in order to address nasal obstruction.

The stent is deployed following repositioning of the cartilages, holds the tissue in place during the healing process, and is removed in approximately 10 days.

The ALAR Nasal Valve Stent is contraindicated for patients with insufficient or over-resected alar cartilages from prior surgery or non-operated, weak or malpositioned lateral crura.



15000

Product	Description	Qty
15000	20 x 16 mm Fluoplastic Stent	2 per box

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications.

For further information, please call Medtronic ENT at 800.874.5797 or consult Medtronic's website at www.medtronic.com/ent.

Medtronic

Medtronic plc ENT

6743 Southpoint Drive N
Jacksonville, FL 32216
USA
Toll free: (800) 874-5797
Telephone: (904) 296-9600
Fax: (800) 678-3995

International Telephone Numbers

Europe HQ Switzerland 41-21-802-7000
Latin America HQ 305-500-9328
Adriatic Region 385-1-488-1120
Argentina 54-11-4898-5700
Australia 1-800-668-670
Baltic Region 37-1-67560226
Belgium 32-2456-09-09
Brazil 55-11-2182-9200
Canada 1-800-268-5346
Chile 56-2-2655-5110
China 86-10-5869-8989
Colombia 57-1-742-7300
Czech Republic 420-233-059-111
France 33-470-679-800
Germany 49-2159-8149-353
Greece 30-210-67-79-099
Hong Kong 852-2919-1300
Hungary 36-188-90600

India 91-22-33074700
Israel 972-9-972-4400
Italy 39-02-24137-324
Japan 81-3-6774-4611
Korea 82-2-3404-3600
Lebanon 961-1-370-670
Luxembourg 32-2456-09-09
Malaysia 60-37-953-4800
Mexico 52-55-11-02-90-30
Netherlands 31-45-566-8800
Poland 48-22-4656900
Russian Fed. 7-495-580-73-77
Singapore 65-6436-5000
South Africa 27-11-260-9300
Spain 34-91-625-05-40
Taiwan 886-2-2183-6000
Thailand 662-232-7400
UK 44-1923-205-166