





















## لیست هندپیس های لیگاشور برای عمل های باز

					
<b>MODEL</b>	<b>LF4418</b>	<b>LF2019</b>	<b>BZ4212</b>	<b>LF1923</b>	<b>LS1020</b>
<b>Description</b>	<b>LigaSure Impact™ Open Instrument</b>	<b>LigaSure™ Exact™ Open Instrument</b>	<b>LigaSure™ Bizact™ Tonsillectomy</b>	<b>LigaSure™ Maryland Jaw 23 cm Open Instrument (Nano Coated)</b>	<b>LigaSure Atlas™ 20 cm Open Instrument</b>
<b>Tip</b>					
<b>Length</b>	18 cm (7.1 in)	18.8 CM	12 CM	23 CM	20 CM
<b>Shaft Diameter</b>	13.5 mm (oval)	-	-	5 MM	10 MM
<b>Shaft Rotation</b>	180 Degrees	-	-	360 DEGREES	359 DEGREES
<b>Jaw Design</b>	Curved, 14 Degrees	-	-	Maryland jaw	Straight, Bilateral, Blunt
<b>Seal Length</b>	36 mm	16.5 MM	14 MM	20 MM	22 MM
<b>Cut Length</b>	34 mm	14.7 MM	12 MM	18 MM	20 MM
<b>Procedural applications</b>	<ul style="list-style-type: none"> <li>❖ Urology</li> <li>❖ Colorectal</li> <li>❖ General surgery</li> <li>❖ Gynecology</li> </ul>	<ul style="list-style-type: none"> <li>❖ Ear, Nose and Throat (ENT)</li> <li>❖ General</li> <li>❖ Plastic/Reconstructive</li> <li>❖ Urologic</li> <li>❖ Thoracic</li> </ul>	<ul style="list-style-type: none"> <li>❖ Tonsillectomy</li> </ul>	<ul style="list-style-type: none"> <li>❖ Adhesiolysis</li> <li>❖ Abdominal hysterectomy</li> <li>❖ Gastric bypass</li> <li>❖ Colon resection</li> <li>❖ Cystectomy</li> <li>❖ Radical prostatectomy</li> <li>❖ Gastrectomy</li> <li>❖ Splenectomy</li> <li>❖ Salpingo-oophorectomy</li> <li>❖ Nephrectomy</li> </ul>	<ul style="list-style-type: none"> <li>❖ Adhesiolysis</li> <li>❖ Abdominal hysterectomy</li> <li>❖ Gastric bypass</li> <li>❖ Colon resection</li> <li>❖ Cystectomy</li> <li>❖ Radical prostatectomy</li> <li>❖ Gastrectomy</li> <li>❖ Splenectomy</li> <li>❖ Salpingo-oophorectomy</li> <li>❖ Nephrectomy</li> </ul>

## لیست هندپیس های لیگاشور برای عمل های لاپاروسکوپی

					
<b>Model</b>	<b>LF5637</b>	<b>LF1837</b>	<b>LF1937</b>	<b>LS1037</b>	<b>LF1944</b>
<b>Description</b>	<b>LigaSure™ L-HOOK™ 5 Instrument in 1 (Nano Coated)</b>	<b>LigaSure™ 5 mm Blunt Tip 37 cm Laparoscopic Instrument (Nano Coated)</b>	<b>LigaSure™ Maryland Jaw 37 cm Laparoscopic Instrument (Nano Coated)</b>	<b>LigaSure Atlas™ 37 cm Hand Switching Laparoscopic Instrument</b>	<b>LigaSure™ Maryland Jaw 44 cm Laparoscopic Instrument (Nano Coated)</b>
<b>Tip</b>					
<b>Length</b>	37cm	37cm	37cm	37 CM	44 CM
<b>Shaft Diameter</b>	5mm	5mm	5 MM	10 MM	5 MM
<b>Shaft Rotation</b>	350 Degrees	180 degrees	350 DEGREES	359 DEGREES	360 DEGREES
<b>Jaw Design</b>	Maryland LigaSure with L-Hook Electro-surgery	Straight, Blunt	Maryland jaw	Straight, Bilateral, Blunt	Maryland jaw
<b>Seal Length</b>	19.5 mm	19.5 mm	20 MM	22 MM	20 MM
<b>Cut Length</b>	17.8 mm	17.8 mm	18 MM	20 MM	18 MM
<b>Procedural applications</b>	<ul style="list-style-type: none"> <li>❖ Adhesiolysis</li> <li>❖ Adrenalectomy</li> <li>❖ Colectomy</li> <li>❖ Gastrectomy</li> <li>❖ Gastric bypass</li> <li>❖ Laparoscopic hysterectomy</li> <li>❖ Nephrectomy</li> <li>❖ Nissen fundoplication</li> <li>❖ Oophorectomy</li> <li>❖ Splenectomy</li> </ul>	<ul style="list-style-type: none"> <li>❖ Adhesiolysis</li> <li>❖ Adrenalectomy</li> <li>❖ Colectomy</li> <li>❖ Laparoscopic hysterectomy</li> <li>❖ Nephrectomy</li> <li>❖ Oophorectomy</li> <li>❖ Roux-en-Y gastric bypass</li> <li>❖ Sleeve gastrectomy</li> <li>❖ Splenectomy</li> <li>❖ Whipple</li> </ul>	<ul style="list-style-type: none"> <li>❖ Adhesiolysis</li> <li>❖ Adrenalectomy</li> <li>❖ Colectomy</li> <li>❖ Gastrectomy</li> <li>❖ Gastric bypass</li> <li>❖ Laparoscopic hysterectomy</li> <li>❖ Nephrectomy</li> <li>❖ Nissen fundoplication</li> <li>❖ Oophorectomy</li> <li>❖ Splenectomy</li> </ul>	<ul style="list-style-type: none"> <li>❖ Adrenalectomy</li> <li>❖ Axillary dissection</li> <li>❖ Gastrectomy</li> <li>❖ Colectomy</li> <li>❖ Nephrectomy</li> <li>❖ Liver resection</li> <li>❖ Splenectomy</li> <li>❖ Whipple</li> </ul>	<ul style="list-style-type: none"> <li>❖ Adhesiolysis</li> <li>❖ Adrenalectomy</li> <li>❖ Colectomy</li> <li>❖ Laparoscopic hysterectomy</li> <li>❖ Nephrectomy</li> <li>❖ Oophorectomy</li> <li>❖ Roux-en-Y gastric bypass</li> <li>❖ Sleeve gastrectomy</li> <li>❖ Splenectomy</li> <li>❖ Whipple</li> </ul>