




Ova-Stiff™ Ovum Aspiration Needles—B Bevel

Used for ovum aspiration.

- Improved Chiba-type bevel (B bevel) grind creates the sharpest needle available to reduce patient discomfort.
- Larger ergonomic handle for greater comfort and improved control during use.
- Extra-stiff and extra-smooth needle cannula enables precise placement of the needle tip in follicle.
- EchoTip® enhances visualization of needle tip when used with ultrasonic imaging equipment.




Global Product Number	Order Number	gage	Length cm	Aspiration Line cm
G30430	K-OSN-1630-B-60	16	30	60
G30432	K-OSN-1635-B-60	16	35	60
G29617	K-OSN-1635-B-90	16	35	90
G29618	K-OSN-1730-B-60	17	30	60
G29625	K-OSN-1735-B-60	17	35	60
G30460	K-OSN-1730-B-90	17	30	90
G29626	K-OSN-1735-B-90	17	35	90



Immature Ovum Aspiration Needle

Used for laparoscopic or ultrasound-guided transvaginal aspiration of oocytes from small ovarian follicles. Smaller needle gages are also suitable for immature ovum aspiration.

- EchoTip® extends around the bevel, improving visibility and placement of the needle.
- 17 gage, 30 cm guide needle used to puncture vaginal wall, adding rigidity to the 20 gage, 35 cm aspiration needle.
- New stopper design accepts the tip of a Luer slip syringe, facilitating follicle flushing.



Global Product Number	Order Number	Guide Needle gage	Guide Needle Length cm	Aspiration Needle gage	Aspiration Needle Length cm	Aspiration Tubing Length cm
G38609	K-IOPS-2035-1730	17	30	20	35	60