



Giving Patients Their Voice Back.

VOX Implants™

Injectable PDMS Elastomer for Vocal Cord Rehabilitation



VOX Implants™

■ Injectable PDMS Elastomer for Vocal Cord Rehabilitation

VOX Implantation procedure is safe, simple and effective method of restoring natural speech

By producing medial displacement to allow apposition by the normal vocal cord , VOX Implantation may be an ideal primary treatment for both unilateral laryngeal paralysis or correction of deglutition difficulties after laryngeal reconstructive surgery. Studies have demonstrated that VOX Implants are biocompatible and non-migratory. VOX procedure creates a soft, permanent mucosal elevation.

- Long-lasting results
- Minimally invasive surgery
- Non-Migratory
- Biocompatible and non-absorbable
- Top-quality, reusable application instrument

VOX Implants consist of latex-free, non-pyrogenic, solid, irregularly-textured, medical-grade silicone elastomer implants, suspended in a hydrogel carrier, and packaged in a clear polycarbonate syringe.

Ordering code	Description
VOX-1.0	VOX Implant
VMN-720	VOX Malleable Needle
BP-AD	VOX Administration Instrument

www.medtronic.eu

Europe

Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
www.medtronic.eu
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00